



STEP 1

Release the tie straps



STEP 2

Straighten stretcher collapsible support rods



STEP 3

Hold both handles and rotate 90° outwards.

Dimensions Extended
2290mm X 584mm X 171mm



STEP 4

While holding handles, push against cross rod and lock into place.

Turn stretcher around and lock second rod into place.



Copyright Golden Season Pte Ltd. All Rights Reserved.



Golden Season Pte Ltd | 31 Playfair Road, Singapore 367993
T 65.6863 2228 | F 65.6363 2478
Email: info@goldenseason.com.sg
Website: www.goldenseason.com.sg | www.gseason.com

Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2