



GS3 Ambulance Stretcher



The GS3 Ambulance Stretcher is designed for patient evacuation and transportation by ambulance services and hospitals.

The legs of the stretcher fold down to allow entry into an ambulance stretcher, and thereafter locks into position with the bracket (provided with product).

The GS3 is equipped with a gas spring, allowing for varied resting angles from 0 to 90 degrees. The adjustment can be done even when the patient is on the stretcher without the need to move the patient and risk further trauma.

It is also equipped with collapsible side rails to facilitate placement of patient onto the stretcher. The product comes complete with locking bracket, resting cushions and IV stand.

Golden Season Pte Ltd



Key features at glance:

- Gas spring operated backrest
- Collapsible side rails
- Folding Legs
- Rugged and durable
- Easy to operate

GS3 Ambulance Stretcher – Specifications



Specifications:	
Length:	190 cm
Width:	54 cm
Height:	96 cm
Low position height:	26 cm
Resting angle:	0° to 90°
Maximum Recommended Load:	136 kg
Weight:	44.5 kg

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478
Email: info@goldenseason.com.sg Website: www.goldenseason.com.sg | www.gseason.com Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2