



GS5 Ambulance Stretcher



The GS5 Ambulance Stretcher is designed for city conditions where patients are commonly found on high rise apartments.

The main feature of the product is it's ability to fold down into a wheelchair, allowing for access in small elevators. The conversion from Ambulance stretcher to wheelchair can be done even when the patient is on the stretcher.

The GS5 is also equipped with a gas spring to allow varied resting angles for the patient. It can be operated even when the patient is on the stretcher, preventing unnecessary trauma associated with movement with injuries.

The legs of the stretcher fold down to allow entry into ambulances and locks with the bracket provided.

The product comes complete with IV drip stand, locking bracket, resting cushions, lockable wheels.

Key features at glance:

- Able to convert to wheelchair
- Rugged and Durable
- Gas spring operated backrest
- Collapsible side rails
- Easy to operate



GS5 Ambulance Stretcher - Specifications



Specifications	
Length:	190cm
Width:	54cm
Height:	89cm
Low position height:	25cm
Resting angle:	0° to 80°
Max. Load:	159kg

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478
Email: info@goldenseason.com.sg Website: www.goldenseason.com.sg | www.gseason.com Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2