



MedEvac4 Folding EMS Stretcher

Manufactured with



Specifications	
Dimension (Extended)	Approximately 2150mm x 565mm x 100mm
Dimension (Folded)	Approximately 547mm x 216mm x 190mm
Weight	Approximately 6.5kg
Frame material	Altrunium® with powder coating finish
Canvas type	Heavy Duty Laminated PVC, waterproof
Canvas specifications	Weight: 240gsm Warp: Minimum 2400N / 50mm Weft: Minimum 1400N / 50mm
Patient restraints	Inbuilt to canvas with buckle
IV pole attachment points	6 slots for our Tactical Adjustable IV pole
Handle type	Ergonomically grooved to human hand, black
Payload	570kgs without material fatigue
Static load	Over 750kgs
Midpoint deflection with 95kg	load Less than 2cm
Dynamic presses	Over 500,000 repetitions with 250kg presses
Compatible accessories	MedEvac4 Litter Stands, IV pole, MedEvac4 Backpack, MedEvac4 Backrest, MedEvac4 Wheelcart
The MedEvac4 Stretchers are compliant with Military Standards MIL-L-16446A	

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd