



Scoop Stretcher



The Golden Season Scoop Stretcher is designed for transportation of patients suffering from back injuries without adding additional trauma. Standard evacuation platforms contact the patient's spine adding pain and discomfort during transportation.

Scooping up the patient

The scoop stretcher is able to effectively scoop up the patient from laying position.

Foldable design

Our scoop stretchers are foldable and compact for storage. It is easy to transport as well.

Extremely Durable

Golden Season's scoop stretchers utilize rugged joints and an adamant construction to withstand abuse and rough handling.

Specifications	
Dimension when expanded	170cm x 43cm x 9cm
Dimension when folded	85cm x 43cm x 11cm
Weight:	Approximately 10 kg
Accessories:	2 / 3 patient restraints provided

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd