



MedEvac2 Tactical Stretcher

ALRUNIM[®]



FEATURE

The MedEvac 2 combat stretcher is the cost-effective casualty platform for rescuers faced with the challenges of natural disasters and mass casualty evacuation. Bifold, light weight with extremely high tensile strength. Able to deploy quickly and efficiently.

Super Compact

The MedEvac2 is a lightweight folding stretcher that folds into two segments. Only 42.5" (1080cm) when folded. Ideal to be deployed in field hospitals and evacuation centres

Lightweight

The MedEvac2 weights only 5.5kg (12.2lbs). It opens fast and easily (within 6 seconds). Just snap it open and its patented joint configuration automatically locks in place. No need to waste precious seconds fiddling with locking mechanisms as with other folding stretchers. Once it opens, it stays open - no chance of accidents - the automatic locking joints guarantee it. Refolding is just as fast - only seven seconds.

Golden Season Pte Ltd



Not just for military

All the characteristics that make the MedEvac2 the ideal field hospital stretcher also make it perfect for a wide range of civilian rescue and evacuation operations. Firefighting, helicopter evacuation, ski patrol, wilderness expeditions, camping, etc.

The MedEvac2 is the best stretcher for the job; lightweight, portable, easy to use and reliable.

The MedEvac2 tactical stretcher is the best stretcher for the job; lightweight, portable, easy to use and reliable.

Rugged and reliable simplicity

The MedEvac2's revolutionary joint configuration eliminates the need for springs and complicated locking mechanisms. With only 2 folding joints on either side, the MedEvac2 opens and closes using natural motions. It doesn't take a genius to use the MedEvac2, only to design it. The poles are tested to take 570Kgs of stress without signs of material fatigue.

Easy transportation and storage

The MedEvac2's small folded volume and light weight means more beds can be stowed in tight places. That's important in evacuation carriers such as helicopters, where space and weight is at a premium.

Built from ALTRUNIUM™



The MedEvac2 Tactical Stretcher is built using ALTRUNIUM® - an extremely stable molecular structure formulated by our team of experienced metallurgist.



MedEvac2 Tactical Stretcher – Specifications

Technical Specifications	
Dimension when Extended:	2140mm X 560mm X 145mm
Dimension when Folded:	1085mm X 195mm X 190mm
Weight:	5.5kgs
Frame Material:	Altrunium® with powder coating finish
Canvas Specifications:	Weight: 240gsm Warp: Minimum 2400N / 50mm Weft: Minimum 1400N / 50mm
Patient Restraints:	Built-in patient restraints with buckle
IV Pole attachment points:	2 slots for our Tactical Adjustable IV Pole
Cross bar support:	Reinforced ring made from Altrunium®
Handle Type:	Grooved to human hand tactical black
Payload:	570kg without material fatigue
Midpoint Deflection with 95kg load:	Less than 2cm
Ground Clearance:	10cm
Static Load:	Over 750kg
Dynamic presses with:	Over 50,000 repetitions with 250kg presses
Compatible Accessories:	MedEvac4 Litter Stands, IV pole, MedEvac4 Backpack, MedEvac4 Backrest, MedEvac4 Wheelcart
Leg crossbar:	Altrunium® crossbar with powder coating finish and support straps to prevent collapsing.

Accessories	
Transfusion Bag Hanger (IV Pole)	Fits into holes provided
Litter Stand	Elevation provides doctors and examiners a better position to examine the patient, as well as perform field surgery
First Aid Station	A Combination of MedEvac2 Tactical Stretcher, Tactical IV Pole, Blanket and First Aid kit It can be ordered in a variety of configurations to suit specific applications such as airports, public buildings, factories etc..
MedEvac2 Stretcher Backrest	Designed for MedEvac2 Tactical Stretcher that conforms to the patient body contour to relieve discomfort

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478
Email: info@goldenseason.com.sg Website: www.goldenseason.com.sg | www.gseason.com Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2



MedEvac2 Tactical Stretcher – Accessories

Altrunium® Reinforced Ring



IV pole Attachments, 2 points per stretcher Compatible with our Tactical IV Pole



Altrunium® Leg Crossbar with Strap support to prevent unintentional collapsing on impact.



The MedEvac4 is manufactured to military specification Mil-L-16446A and Certified by Third Party test centers for the above specifications. Test reports are available.

Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478
Email: info@goldenseason.com.sg Website: www.goldenseason.com.sg | www.gseason.com Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2



MedEvac2 Tactical Stretcher – Accessories



The MedEvac2 Tactical Stretcher can be used with the following compatible products to enhance field capabilities and scope of use:



MedEvac4 Stretcher Backrest



MedEvac4 Tactical Litter Stand



MedEvac4 IV Pole



Stretcher Wheel Cart