

MEDEVAC4 STRETCHER

STANAG 2040 NATO Stretcher

The MedEvac4 Fold STANAG 2040 Stretcher is the NATO version of the reliable and proven MedEvac4 Stretcher Platform, and is manufactured in accordance with NATO Standards [STANAG 2040]. It folds into 4 for compact transport / storage, and yet maintains the extremely high load bearing capability that our stretchers are renown for.

This appliance conforms to the Directive 93/42/CEE “Medical Devices”



Specifications

Dimensions Extended	2290mm X 584mm X 171mm
Dimensions Folded	630mm X 215mm X 220mm
Canvas Specification	Warp: Minimum 2400N / 50mm Weft: Minimum 1400N / 50mm Compliant with MIL-L-16446A Standards
Frame	Aerospace Grade Aluminum
Handles	Retractable, Automatic locking, Ergonomically Designed
Weight	7kg
CE Mark	Yes

Specifications are subject to change without prior notice.

Key Features & Materials

- Aerospace Grade Aluminum
- Automatic Locking Handles
- Fast & Easy to assemble
- Compliant with MIL-L-16446A Standards



Copyright Golden Season Pte Ltd. All Rights Reserved.

