



Stretcher Wheel Cart/Wheeled Litter Carrier

(PATENTED)



Golden Season's Stretcher Wheel Cart/Wheeled Litter Carrier is a lightweight casualty transport system that allows one person to transport a casualty strapped to a stretcher, and also performs as a litter stand in the battlefield or medical stations.

The revolutionary design is foldable, and can also be totally disassembled for compact storage, easy shipping and transportation. The Stretcher Wheel Cart/Wheeled Litter Carrier is compatible with modern combat litters, and is equipped with non puncture wheels. A single user can operate the Stretcher Wheel Cart/Wheeled Litter Carrier.



Technical Information

Deployed dimensions:	L143 x W70 x H92 (cm)
Half Folded dimensions:	L126 x W70 x H50 (cm)

Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478
Email: info@goldenseason.com.sg Website: www.goldenseason.com.sg | www.gseason.com Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2



Fully Folded dimensions:	L107 x W70 x H25 (cm)
Weight:	17.5kgs excluding litter

*Specifications are subject to change without prior notice.



Characteristics

- ✓ Foldable and compactable
- ✓ Able to perform as a litter stand
- ✓ Lightweight – 17.5kgs only
- ✓ Pneumatic non puncture wheels
- ✓ Can be fully disassembled for easy shipping

The Stretcher Wheel Cart/Wheeled Litter Carrier is compatible with modern stretcher systems including the Medevac4, TalonII 90C, Raven 90C, T3 Assault Litters and other stretchers.

Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478
Email: info@goldenseason.com.sg Website: www.goldenseason.com.sg | www.gseason.com Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2