



## MedEvac4 Backpack

MedEvac4 Backpack



The MedEvac4 Backpack is designed to carry 1 X MedEvac4 Tactical Stretcher, as well as medical supplies and first aid equipment used by the combat medic. It is ergonomically designed for fast access with slots for labeling. The back of the stretcher is outfitted with 2 metal strips to prevent the stretcher from knocking against the spine of the combat medic when running.



Please email us at [info@goldenseason.com.sg](mailto:info@goldenseason.com.sg) for any queries relating to the product above.

### Golden Season Pte Ltd

31 Playfair Road Singapore 367993 T 65.6863 2228 F 65.6363 2478

Email: [info@goldenseason.com.sg](mailto:info@goldenseason.com.sg) Website: [www.goldenseason.com.sg](http://www.goldenseason.com.sg) Co. Regn No: 03734/1984E

GST Registration No. M2-0067495-2