



MK II MedEvac4 STANAG 2040 NATO Stretcher



The MedEvac4 Fold STANAG 2040 Stretcher MK II is improved from the original version of the reliable and proven MedEvac4 Stretcher Platform, and is manufactured in accordance with NATO Standards (STANAG 2040). The lightweight Altrunium® frame folds by 4 for compact transport / storage, and maintains the extremely high load bearing capability that our stretchers are renown for. The stretcher has 6 built in IV Pole attachment points to support intravenous drip bags, and a waterproof, fire retardant mesh canvas.

CE This appliance conforms to the Directive 93/42/CEE “Medical Devices”

Specifications:	
Dimensions Extended	Approximately 2286mm X 575mm X 168mm
Dimensions Folded	Approximately 515mm X 230mm X 240mm
Canvas Specification	Waterproof, fire retardant perforated mesh. Easy to disinfect and clean, does not react to bleach. Load tested to 1225kgs loaded evenly over stretcher
Frame	Altrunium® Frame, load tested over 1000kgs loaded evenly
Handles	Retractable, Automatic locking, Ergonomically Designed
Weight	7.85 kgs
CE Mark	Yes

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478
Email: info@goldenseason.com.sg Website: www.goldenseason.com.sg | www.gseason.com Co. Regn No: 03734/1984E
GST Registration No. M2-0067495-2